



Email : sales@laboratoryinstrumentindia.com

Product Name :
Complex Resharpener Of Mill and Drill

Product Code :
LBNY-0005-12300051



Description :

Grinding Wheel: CBN Wheel For HSS Drill Bit/SD Wheel for Carbide Drill Bit.
Standard equipment: Twelve collets, Nineteen collets.
Complex re-sharpener of end mill & drill, one machine with two function, easy to use, economical and practical, with no skill requirement to grind.

Technical Specification :

Complex Resharpener Of Mill and Drill



Laboratory Instrument India