



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Deluxe Model of Shaping Machine

**Product Code :**  
LBNY-0005-12300053



**Description :**

**Features:**

Specially designed shaper machine for increased-production, precision Performance and Deep cut.  
Selected special alloy steel used for screws and other steel parts.  
Naturally Seasoned Casting.  
Quick Return Motion.  
Heavily cross ribbed, Robust and sturdy construction ensure free working.  
Fully Hand Scrapped.

**Technical Specification :**

Deluxe Model of Shaping Machine



**Laboratory Instrument India**