



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
AC Power Supply

**Product Code :**  
LBNY-0002-610001



#### **Description :**

AC Power Supply, AC Power Supply Equipments, AC Power Supply Tools, AC Power Supply Tool Kits, AC Power Supply Manufacturers, AC Power Supply Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia

#### **Technical Specification :**

The loads are designed for testing UPS devices, inverters, AC sources, and other components such as switches and fuses that are disconnected from the mains.

The loads can simulate load conditions with a high crest factor and different power factors with real-time compensation, even with a distorted voltage waveform.

This special function enables a simulation of the real world and prevents overloading, which provides reliable and accurate test results.

Digital display for output voltage and current

Output voltage and current continuously adjustable

Over current protection

Isolative safety

Input voltage: 230 (220) VAC

---

Output voltage: 0~260 VAC  
Protection: over current.



**Laboratory Instrument India**