



Email : sales@laboratoryinstrumentindia.com

Product Name :
Workshop Milling Machine Parts

Product Code :
LBNY-0004-760002



Description :

Workshop Milling Machine Parts

Technical Specification :

Milling machine
Magnetic Contactor: 40A,380V
Thermal Overload: 16A
Three Phase: Mcb, 32a 32A
Single Phase Mcb: 10,16, 6,
Pushbutton: RED, GREEN
Auxiliary Contact: 10A,
Din Rail: 1.5
Phase Indicator: 380V
Single Core Cable, Red: 2.5mm²
Industrial extension cable: 2.5mm, 3-core
250V, 100metres with (4pin), female outlet

Workshop Lathe Machine parts, Workshop Lathe Machine parts Equipments, Workshop Lathe Machine parts Tools, Workshop Lathe Machine parts Tool Kits, Workshop Lathe Machine parts Manufacturers, Workshop Lathe Machine parts Suppliers from India, China, Kenya



Laboratory Instrument India