



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Workshop Milling Machine Parts

**Product Code :**  
LBNY-0004-760002



**Description :**

Workshop Milling Machine Parts

**Technical Specification :**

Milling machine  
Magnetic Contactor: 40A,380V  
Thermal Overload: 16A  
Three Phase: Mcb, 32a 32A  
Single Phase Mcb: 10,16, 6,  
Pushbutton: RED, GREEN  
Auxiliary Contact: 10A,  
Din Rail: 1.5  
Phase Indicator: 380V  
Single Core Cable, Red: 2.5mm<sup>2</sup>  
Industrial extension cable: 2.5mm, 3-core  
250V, 100metres with (4pin), female outlet

Workshop Lathe Machine parts, Workshop Lathe Machine parts Equipments, Workshop Lathe Machine parts Tools, Workshop Lathe Machine parts Tool Kits, Workshop Lathe Machine parts Manufacturers, Workshop Lathe Machine parts Suppliers from India, China, Kenya



**Laboratory Instrument India**