



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Ball End Mill Resharpener

**Product Code :**  
LBNY-0005-12300054



**Description :**

With Taiwan diamond grinding wheel, only one piece can complete all procedure.  
It can be equipped directly with an accurate angle and long service life.  
It is the special grinder for ball end mill.  
Grinding is accurate and rapid, easy operation, with no skill requirement to grind.

**Technical Specification :**

Grinding Range: 6-20  
Motor: 220V/(250W)  
Speed: 4800rpm  
Grinding Wheel: CBIM(HSS)orSD  
Dimension: 60 x 29 x 37cm  
Weight: 36.5kg



**Laboratory Instrument India**