



Email : sales@laboratoryinstrumentindia.com

Product Name :
Ball End Mill Resharpener

Product Code :
LBNY-0005-12300054



Description :

With Taiwan diamond grinding wheel, only one piece can complete all procedure.
It can be equipped directly with an accurate angle and long service life.
It is the special grinder for ball end mill.
Grinding is accurate and rapid, easy operation, with no skill requirement to grind.

Technical Specification :

Grinding Range: 6-20
Motor: 220V/(250W)
Speed: 4800rpm
Grinding Wheel: CBIM(HSS)orSD
Dimension: 60 x 29 x 37cm
Weight: 36.5kg



Laboratory Instrument India