



Email : sales@laboratoryinstrumentindia.com

Product Name :
Drill Bit Resharpener

Product Code :
LBNY-0005-12300056



Description :

Grinding Wheel: CBN Wheel For HSS Drill Bit / SD Wheel for Carbide Drill Bit

Portable drill grinder. With Taiwan Grinding wheel, it can be equipped directly with an accurate angle and long service life.

It is not so easy to align the drill less than 3mm.

Drill Bit Resharpener can sharpen the drill bit front cutting lip, lip relief angle, you also can control the center spot at random instead of center drill, escape of chips easily, drill with a light heart.

Technical Specification :

Drill Bit Resharpener



Laboratory Instrument India