



Email : sales@laboratoryinstrumentindia.com

Product Name :
Drill Bit Resharpener

Product Code :
LBNY-0005-12300058



www.laboratoryinstrumentindia.com

Description :

Grinding with travel block control ensures the bit high accuracy, symmetry, and not only the penetrating power, but also longer bit life. Drill Bit Re-Sharpener drill bit grinder, with the camera to set the drill grinding and do the measure after grinding. It can grind drill bit with flat relief angle and various point thinning.

Technical Specification :

Range: 1-6
Point Angle: 90-150 deg.
Power: 250W
Speed: 4800rpm
Dimension: 62x34x33cm
Weight: 28kg/30kg



Laboratory Instrument India