



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Standard Model of Slotting Machine

**Product Code :**  
LBNY-0005-12300059



**Description :**

Slotting Machine can produce slotting strokes from 10 to 150mm. This equipment is made of rigid construction and it creates no vibration while machining. It is made of fine grained and seasoned casting that makes the instrument work for long term application. The Ram adjustment is 150mm and the column to Ram distance is 210mm similarly table to Ram distance is 310 mm. The longitudinal and cross movement for ranges about 200mm and 110 mm respectively.

**Technical Specification :**

ADJUSTABLE Stroke (mm): 10 to 150 10 to 250  
Longitudinal Movement (mm): 200 230  
Cross Movement (mm): 110 230  
Speed adjustment (Speed): 3 3  
Ram Adjustment (mm): 150 250  
Motor (HP): 1 1.5  
Column to RAM Dist.(mm): 210 330  
Table to RAM Dist. (mm): 310 410  
Table size (mm): 280x335 280x335  
Machine Dimension: 1050x550x1175 1150x600x1375  
Machine weight (kg): 375 450

---

# Laboratory Instrument India

