



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Screw Tap Resharpener

**Product Code :**  
LBNY-0005-12300062



**Description :**

Range: M5-M20 M16-M42  
Point Angle: 5 -30 deg. 5 -30 deg.  
Power: 220V 180W 220V 250W  
Motor: 4800rpm 4800rpm  
Grinding wheel: CBN(HSS)orSD  
Dimension: 36x27x30.5cm 40x30x31cm  
Weight: 16kg 27kg

**Technical Specification :**

Screw Tap Resharpener



**Laboratory Instrument India**