



Email : sales@laboratoryinstrumentindia.com

Product Name :
Capstan Lathe

Product Code :
LBNY-0005-12300063



Description :

The turret head is flexible for making sliding movement over the ram. It stays suitable equipment for industries which require identical kind of tooling frequently. The Capstan lathe belongs to the light duty kind of lathes and it is featured with turret head on the ram. This equipment can also be applied to obtain different cuts by rotating the turret head on the block wherein it becomes quite flexible and fast to work with. The ram is then fixed on the saddle which will not move at the time of tooling. It suits perfect for performing machining over short work materials. Above all, when it comes to user friendly features in cutting, Capstan apparatus from Machine Tools is one of the best options.

Technical Specification :

Capstan Lathe



Laboratory Instrument India