




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Standard Model Of Shaping Machine

**Product Code :**  
LBNY-0005-12300064



**Description :**

**Features:**

Specially designed for increased-production, precision Performance and Deep Cut.  
Fully Hand Scrapped  
Naturally Seasoned Casting.  
Selected special alloy steel used for screws and other steel parts.  
Heavily cross ribbed, Robust and sturdy construction ensures free working.  
Quick Return Motion.

**Technical Specification :**

Standard Model Of Shaping Machine



**Laboratory Instrument India**