




Email : sales@laboratoryinstrumentindia.com

Product Name :
Standard Model Of Shaping Machine

Product Code :
LBNY-0005-12300064



Description :

Features:

Specially designed for increased-production, precision Performance and Deep Cut.
Fully Hand Scrapped
Naturally Seasoned Casting.
Selected special alloy steel used for screws and other steel parts.
Heavily cross ribbed, Robust and sturdy construction ensures free working.
Quick Return Motion.

Technical Specification :

Standard Model Of Shaping Machine



Laboratory Instrument India