



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Tool and Cutter Grinding Machine

**Product Code :**  
LBNY-0005-12300065



**Description :**

An adjustable stopper limits the longitudinal table reciprocation movement.  
The provision of tilting head simplifies the cutting face angle quickly and accurately  
Work head holding fixture is fully adjustable with swivel base.  
Almost all cylindrical side & face, involute gear, Concave as well as convex cutters,  
Drill pints, tap lead & Radius grinding jobs can also be done without problem with extra accessories.  
Bevel lead angle reamers & Cutting face of tap flute can be resharpened with standard accessories.  
The tool runs at mono speed of 3600 RPM and the drive from the motor to spindle is by means of v-belt.

**Technical Specification :**

Tool and Cutter Grinding Machine



**Laboratory Instrument India**