



Email : sales@laboratoryinstrumentindia.com

Product Name :
Double Column Semi Automatic Horizontal Bandsaw
Machine

Product Code :
LBNY-0005-12300071



Description :

Silent Features:

- Automatic Height Adjustment for Saw Frame
- Variable Pressure Control for Cutting Feed
- Stress Relieved Fabricated Steel Structure
- Heavy Duty Worm Gear Box
- Moveable Blade Support Guide
- Manual Blade Tensioning Arrangement
- Hard Chrome Plated Ground Column
- Auto Cutting Stop
- Sawing mode changes, flexible, manual & automatic Combined
- PLC controller complete automatic sawing process
- Push Button Type Operating Panel Board
- Material Feeding Manually.

Technical Specification :

Double Column Semi Automatic Horizontal Bandsaw Machine

Laboratory Instrument India

