



Email : sales@laboratoryinstrumentindia.com

Product Name :
Drilling And Milling Machine Auto Feed

Product Code :
LBNY-0005-12300075



Description :

Features:

Spindle speed step-change cyclically, spindle auto feeding and cyclic change of feeding speed
Precision micro feed
Positive spindle lock
Distinct auto device for releasing tools, operate very easily
Gear driven, lower noise
Milling, drilling, and tapping
Super high column
Head swivels Horizontally.

Technical Specification :

Drilling And Milling Machine Auto Feed



Laboratory Instrument India