



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Drilling And Milling Machine Auto Feed

**Product Code :**  
LBNY-0005-12300075



**Description :**

**Features:**

Spindle speed step-change cyclically, spindle auto feeding and cyclic change of feeding speed  
Precision micro feed  
Positive spindle lock  
Distinct auto device for releasing tools, operate very easily  
Gear driven, lower noise  
Milling, drilling, and tapping  
Super high column  
Head swivels Horizontally.

**Technical Specification :**

Drilling And Milling Machine Auto Feed



**Laboratory Instrument India**