



Email : sales@laboratoryinstrumentindia.com

Product Name :
Column Type V Belt Drive Fine Feed

Product Code :
LBNY-0005-12300078



Description :

Features:

Adjustable gibs on table precision
Positive spindle lock
Belt driven and round column
Strong rigidity power cutting and precise positioning.
Micro feed precision
12 steps speed
Milling, drilling, tapping, boring and roaming
Headstock swivels $\pm 360^\circ$ horizontally.

Technical Specification :

Column Type V Belt Drive Fine Feed



Laboratory Instrument India