



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Column Type V Belt Drive Fine Feed

**Product Code :**  
LBNY-0005-12300078



**Description :**

**Features:**

Adjustable gibs on table precision  
Positive spindle lock  
Belt driven and round column  
Strong rigidity power cutting and precise positioning.  
Micro feed precision  
12 steps speed  
Milling, drilling, tapping, boring and roaming  
Headstock swivels  $\pm 360^\circ$  horizontally.

**Technical Specification :**

Column Type V Belt Drive Fine Feed



**Laboratory Instrument India**