



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Column Type Geared Drive Fine Feed

**Product Code :**  
LBNY-0005-12300079



**Description :**

**Features:**

Micro feed precision  
Adjustable gibs on table precision  
Belt driven and round column  
Strong rigidity power cutting and precise positioning.  
Milling, drilling, tapping, boring and roaming.

**Technical Specification :**

Column Type Geared Drive Fine Feed



**Laboratory Instrument India**