



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Pitch Controlled Tapping Machine

**Product Code :**  
LBNY-0005-12300080



**Description :**

**Features:**

All moving parts Finished to close Tolerance  
Spindle Ram is made out of Alloy Steel.  
Column, Arm, Work Head And Gears are made out of Closed Grain.  
Parts Finished on imported Grinding Machine  
Parts Properly Balanced  
Prior Assembling of Main Gear And Pully.

**Technical Specification :**

Pitch Controlled Tapping Machine



**Laboratory Instrument India**