




Email : sales@laboratoryinstrumentindia.com

Product Name :
Hydraulic Drilling Head

Product Code :
LBNY-0005-12300082



Description :

Max. Drilling Capacity: 10 16 25 35
Spindle Travel: 100 150 150 150
Spindle rpm / 3 S: 1000.2000.3000 1000.1500.2000 Optional Optional
Motor HP: 1 2 2 3
Pneumatic Power Kg/cm²: 20 20 20 2

Technical Specification :

Hydraulic Drilling Head



Laboratory Instrument India