




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Hydraulic Drilling Head

**Product Code :**  
LBNY-0005-12300082



**Description :**

Max. Drilling Capacity: 10 16 25 35  
Spindle Travel: 100 150 150 150  
Spindle rpm / 3 S: 1000.2000.3000 1000.1500.2000 Optional Optional  
Motor HP: 1 2 2 3  
Pneumatic Power Kg/cm<sup>2</sup>: 20 20 20 2

**Technical Specification :**

Hydraulic Drilling Head



**Laboratory Instrument India**