



Email : sales@laboratoryinstrumentindia.com

Product Name :
Radial Drilling Machine R40G

Product Code :
LBNY-0005-12300086



Description :

Column, Arm, Work Head And Gears are made out of Closed Grain.
Spindle Ram is made out of Alloy Steel.
Parts Finished on imported Grinding Machine
All moving parts Finished to close Tolerance
Parts Properly Balanced
Prior Assembling of Main Gear And Pully.

Technical Specification :

Radial Drilling Machine R40G



Laboratory Instrument India