



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Radial Drilling Machine 40mm

**Product Code :**  
LBNY-0005-12300087



**Description :**

The diameter of the pillar is 170mm so that it can offer the perfect support while the machine is in the operational mode. The drilling capacity of this drill machine is 40mm which is perfect for different types of industries.

**Technical Specification :**

Radial Drilling Machine 40mm



**Laboratory Instrument India**