



Email : sales@laboratoryinstrumentindia.com

Product Name :
All Geared Auto Feed Radial Drill Machine

Product Code :
LBNY-0005-12300088



Description :

Features:

Column, Arm, Work Head And Gears are made out of Closed Grain.
Spindle Ram is made out of Alloy Steel.
Parts Properly Balanced
Prior Assembling of Main Gear And Pully
Parts Finished on imported Grinding Machine
All moving parts Finished to close Tolerance.

Technical Specification :

All Geared Auto Feed Radial Drill Machine



Laboratory Instrument India