



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Radial Drill Machine 25mm

**Product Code :**  
LBNY-0005-12300089



**Description :**

**Features:**

All moving parts Finished to close Tolerance  
Column, Arm, Work Head And Gears are made out of Closed Grain.  
Spindle Ram is made out of Alloy Steel.  
Parts Finished on imported Grinding Machine  
Parts Properly Balanced  
Prior Assembling of Main Gear And Pulley.

**Technical Specification :**

Radial Drill Machine 25mm



**Laboratory Instrument India**