



Email : sales@laboratoryinstrumentindia.com

Product Name :
Radial Drill Machine 25mm

Product Code :
LBNY-0005-12300089



Description :

Features:

All moving parts Finished to close Tolerance
Column, Arm, Work Head And Gears are made out of Closed Grain.
Spindle Ram is made out of Alloy Steel.
Parts Finished on imported Grinding Machine
Parts Properly Balanced
Prior Assembling of Main Gear And Pulley.

Technical Specification :

Radial Drill Machine 25mm



Laboratory Instrument India