



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
High Precision Pillar Drilling Machine

**Product Code :**  
LBNY-0005-12300090



**Description :**

**Features:**

Main Spindle And Pinon are Made out of Forged Steel.  
Spindle Ram is made out of Steel.  
Parts Properly Balanced  
All Moving Parts Finished to Close Tolerance.  
Parts finished on imported Grinding M/c.  
Prior Assembling of Main Gear And Pulley.

**Technical Specification :**

High Precision Pillar Drilling Machine



**Laboratory Instrument India**