



Email : sales@laboratoryinstrumentindia.com

Product Name :
Geared Radial Drilling Machine

Product Code :
LBNY-0005-1230091



Description :

Features:

- Used for boring to create holes in different work materials.
- Economic in pricing when compared to other boring instruments.
- Can drill maximum of 38mm.
- Creates less noise and less vibration.
- Accuracy in making hole.
- Provides consistent boring.
- Durable and reliable tool.
- No compromise made in quality of machine components.
- Requires less or low maintenance but provides excellent machining performance.

Technical Specification :

Geared Radial Drilling Machine



Laboratory Instrument India

