



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Radial Drill Machine

**Product Code :**  
LBNY-0005-1230092



**Description :**

**Features:**

All moving parts Finished to close Tolerance  
Spindle Ram is made out of Alloy Steel.  
Column, Arm, Work Head And Gears are made out of Closed Grain  
Parts Finished on imported Grinding Machine  
Main Spindle And Pinon are Made out of Steel.  
Parts Properly Balanced  
Prior Assembling of Main Gear And Pulley.

**Technical Specification :**

Radial Drill Machine



**Laboratory Instrument India**