



Email : sales@laboratoryinstrumentindia.com

Product Name :
Radial Drill Machine

Product Code :
LBNY-0005-1230092



Description :

Features:

All moving parts Finished to close Tolerance
Spindle Ram is made out of Alloy Steel.
Column, Arm, Work Head And Gears are made out of Closed Grain
Parts Finished on imported Grinding Machine
Main Spindle And Pinon are Made out of Steel.
Parts Properly Balanced
Prior Assembling of Main Gear And Pulley.

Technical Specification :

Radial Drill Machine



Laboratory Instrument India