



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Workshop Oxygen Plant Filling Service

**Product Code :**  
LBNY-0004-760006



**Description :**

Workshop Oxygen Plant Filling Service

**Technical Specification :**

Oxygen plant, filling service

Oxygen Cylinder Refilling Plants — with capacities from 8-100s of cylinders per day. Complete plants include a feed air compressor, feed air dryer, oxygen generator, oxygen compressor, and a cylinder filling rack.

The oxygen compressor delivers oxygen at up to 2,200 psig (15,169 kPa or 151.6 barg) to a high pressure manifold capable of filling up to 10 cylinders at a time.

These cylinder refilling plants operate automatically and generate nominal 93% oxygen. For special applications, an optional high purity module can be added to the plant, to increase oxygen concentration to 99%  $\pm$ .5%.

Filling oxygen cylinders with oxygen gas

3/8 inches twin welding Hoses

Diameter: 10mm

Length: 50m, W.P 20bar

Workshop Oxygen Plant Filling Service, Workshop Oxygen Plant Filling Service Equipments, Workshop Oxygen Plant Filling Service Tools, Workshop Oxygen Plant Filling Service Tool Kits, Workshop Oxygen Plant Filling Service Manufacturers, Workshop Oxygen Plant Filling Service Suppliers from India, China, Kenya

---

# Laboratory Instrument India

