



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
PLC Analog Learning System

**Product Code :**  
LBNY-0005-12400016



#### **Description :**

The model PLC Analog Learning System includes an analog application station, cable assembly, Control Logix combination analog module, student learning materials for both theory and lab and a teacher's assessment guide. The PLC Analog Learning System adds PLC Troubleshooting Learning System to teach PLC control of analog input and output devices. Applications include temperature control, liquid level control, and variable speed conveyors just to name a few. Students will learn industry relevant skills including how to operate, interface, program and troubleshoot PLC systems using analog I/O. Analog I/O capability enables PLCs to more precisely control processes by providing variable output signals and receiving variable input feedback from sensors.

#### **Features:**

- Self-Contained Learning Platform
- External Analog I/O Control
- Student Reference Guide
- Interchangeable Application Station
- Optional Online eBooks.

#### **Technical Specification :**

PLC Analog Learning System

---

# Laboratory Instrument India

