



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Workshop Cabinet (For Tools Storage)

**Product Code :**  
LBNY-0004-760008



**Description :**

Workshop Cabinet (For Tools Storage)

**Technical Specification :**

**Special Design:** Designed specifically for woodworking tools.

**Slim Fit:** These tool cabinets fit in narrower places, so you can have great tool storage that works with your space.

**Durable Construction:** Constructed from solid maple and maple veneers.

**High Capacity:** Upper interior and doors are lined with pegboard, maximizing hanging space.

**Great Shelving:** Includes one fixed, and two adjustable, 1"-thick plywood shelves.

**Easy Plug-in:** Includes a mounted plug-mold strip with six outlets, switch, breaker, cord, and plug.

**Multipurpose:** A great addition to your maintenance room or woodworking area. Choose from a selection of models to fit different kinds of tools.

**Secure Storage:** With two sets of locking doors, tools stay protected when not in use.

Workshop Cabinet (For Tools Storage), Workshop Cabinet (For Tools Storage) Equipments, Workshop Cabinet (For Tools Storage) Tools, Workshop Cabinet (For Tools Storage) Tool Kits, Workshop Cabinet (For Tools Storage) Manufacturers, Workshop Cabinet (For Tools Storage) Suppliers from India, China, Kenya

---

# Laboratory Instrument India

