




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Board For Electronic Design

**Product Code :**  
LBNY-0005-12400047



**Description :**

It contains all of the necessary elements for the student to carry out the design and wiring of his own electronic circuits, both analogue and digital.

**Technical features:**

- 2 push buttons with normally open and normally closed outputs
- 8 red logic indicators for high levels and 8 green logic indicators for low levels
- 1 Potentiometer - 1 kW
- 8 logical switches
- 1 loudspeaker, 8 W
- 1 Potentiometer - 10 kW
- 2 slide switches with free terminals
- 2 7-segments LCD displays with decoder
- Removable breadboard
- 2mm terminals to connect to the breadboard
- Cables with 2mm banana terminations for connection to the breadboard
- Interface for PC.

**Technical Specification :**

Board For Electronic Design



**Laboratory Instrument India**