



Email : sales@laboratoryinstrumentindia.com

Product Name :
Open Microprocessor Trainer

Product Code :
LBNY-0005-12400049



Description :

The kit consists of a module with power supply, with different interfaces and a set of software application. This 32 bit Microprocessor Trainer can make students understanding and becoming familiar with the programming of the microprocessor and its components.

Technical Features:

- General inputs and outputs with at least 64 ports
- At least 8 Interrupt inputs
- Onboard 256 MB Nand Flash
- Onboard 64MB SDRAM and 2MB Nor Flash
- Keyboard with hexadecimal keys
- Expansion bus for connecting external applications
- Development kit
- 7 segments display
- 10 bit AD converter and 8 bit DA converter.

Technical Specification :

Open Microprocessor Trainer

Laboratory Instrument India

