



Email : sales@laboratoryinstrumentindia.com

Product Name :
DC Power Supply

Product Code :
LBNY-0002-610002



Description :

DC Power Supply

Technical Specification :

DC power supplies with 2 channels, 600W output.
In constant voltage and current operations, the power supply provides high regulation and low ripple & noise.
Output off function
Dual independent adjustable DC outputs.
Independent, series and parallel variety output
Constant voltage and constant current
Low ripple and noise
Four LED displays with 3 digits, selectable for
Two sets
Auto tracking output

DC Power Supply, DC Power Supply Equipments, DC Power Supply Tools, DC Power Supply Tool Kits, DC Power Supply Manufacturers, DC Power Supply Suppliers from India, China, Kenya

Laboratory Instrument India

