



Email : sales@laboratoryinstrumentindia.com

Product Name :
Wood Lathe Machine With Leg

Product Code :
LBNY-0004-770002



Description :

Wood Lathe Machine With Leg

Technical Specification :

Heavy Duty Construction- Heavy-duty, precision ground cast iron bed and cast iron legs ensure stability, durability, strength, and minimal vibration.

Variable Speed Control- There is a speed display panel with ON/OFF button and speed control button for you to operate the lathe machine conveniently.

Large Work Capacity- The wood lathe boasts of 42" (1066 mm) distance between centers, swing over bed 16" (406 mm). High concentricity. Adjust the working distance arbitrarily on the slide rail.

Flexible Operation- The lathe head and the tool rest can rotate 360-degree. You can move tailstock to adjust proper distances for work-pieces. Moving the lock lever to adjust or replace the tool rest.

Wide Application- Perfect for various factories and workshops for processing and polishing various kinds of wood to make wood products..



Laboratory Instrument India

