



Email : sales@laboratoryinstrumentindia.com

Product Name :
Motor Board

Product Code :
LBNY-0005-12400078



Description :

Universal speed control system.
Dual-channel encoder
Built-in four-quadrant amplifier
With integrated four-quadrant display
With variable centrifugal mass

Technical Features:

Main Motor/Generator 12 V, 3000 rpm, 1.2 A, 3.2 Ncm
Load to be connected to the secondary Motor/Generator
Dual optical sensor for speed and direction
Shunt to limit and measure the current
Linear H Bridge to have full motor control.

Technical Specification :

Motor Board



Laboratory Instrument India

