



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Motor Board

**Product Code :**  
LBNY-0005-12400078



**Description :**

Universal speed control system.  
Dual-channel encoder  
Built-in four-quadrant amplifier  
With integrated four-quadrant display  
With variable centrifugal mass

**Technical Features:**

Main Motor/Generator 12 V, 3000 rpm, 1.2 A, 3.2 Ncm  
Load to be connected to the secondary Motor/Generator  
Dual optical sensor for speed and direction  
Shunt to limit and measure the current  
Linear H Bridge to have full motor control.

**Technical Specification :**

Motor Board



**Laboratory Instrument India**

