



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Brick Cutting Machine

**Product Code :**  
LBNY-0004-770004



**Description :**

Brick Cutting Machine

**Technical Specification :**

Brick Cutting Machine

Especially designed copper casting mechanical centrifugal water pump features high pumping head, high flow rate, long service life, and is adequate to cooling cutting blade and remove dust.

The pump is equipped with manual-acting effusion function and convenient water supply system ensures sufficient water supply.

Specifications:

Blade size : 14"/350mm

Block Height :127mm/5"

Table Size (L x W):490mm x 280mm

Voltage : 415V/50Hz

Power: 5.5 – 7.5 hp (Over power - no such machine)

Total Weight: 116 – 130 kg

Dimension (L x W x H) :1200 x 700x1520mm

Accessories 4 pcs 14"/350mm

Diamond Blade



**Laboratory Instrument India**