




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Sequential Logic

**Product Code :**  
LBNY-0005-12400088



#### **Description :**

This board allows the study of the following subjects:  
MSI circuits in TTL logic, shift registers, synchronous and asynchronous counters, astable and monostable multivibrators, seven segment displays.

The board is composed of:

- 4 JK/MS flip flop
- 2 D flip flop
- A binary 4 bit synchronous counter
- A Shift register, 8 bit SI SO
- A BCD asynchronous counter
- A BCD synchronous counter
- 2 BCD rotating switches
- 8 capacitors
- 2 BCD 7 segment decoders
- 2 7 segment displays
- 2 linear potentiometers
- A binary 4 bit asynchronous counter
- 2 Monostable multivibrators.

#### **Technical Specification :**

Sequential Logic



**Laboratory Instrument India**