




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Programmable Logic

**Product Code :**  
LBNY-0005-12400091



**Description :**

This board allows the study of RAM memories and circuits using RAMs for storing logic information and for data transfer to a BUS.

The board is composed of:

Coder for hexadecimal keypad

8 bit analogue to digital converter

8 bit digital to analogue converter

4 + 4 buffers with three state output

4 to 10 line decoder

8 bit latch register (type D)

Hexadecimal keypad

RAM 1024 x 4

Multiturn potentiometer for generating voltages between 0 and 10V.

**Technical Specification :**

Programmable Logic

---

# Laboratory Instrument India

