



Email : sales@laboratoryinstrumentindia.com

Product Name :
Dual Bevel Glide Miter Saw

Product Code :
LBNY-0004-770005



Description :

Dual Bevel Glide Miter Saw

Technical Specification :

The Miter Saw wider cross-cuts, enhanced alignment and smooth glide for the life of the tool. The miter saw's 15 Amp motor delivers no-load speeds of 3,800 RPM for quick, accurate cuts in many materials. For added convenience and versatility, the saw is compatible with a variety of 12 In. blades and it has a clear lower guard to provide unobstructed cut-line visibility. The provides easy adjustments with up front bevel and miter controls and detents at common angles. The saw's all-metal bevel lock lever and range selector are located up front to make bevel adjustments easy.

Size: full size

Style : Miter saw

Material : Stainless Steel

Shape : round

Power Source : corded-electric

Voltage : 120 volts

Item Package

Number of Handles : 1

Handle Material : wood, laminets

Cutting Diameter : 12 inches

Included Components : Blade

Assembled Diameter : 12 inches

Gravity-Rise Miter Saw Stand : yes



Laboratory Instrument India