




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Transistor Multi Vibrators

**Product Code :**  
LBNY-0005-124000104



#### **Description :**

The board is divided into three main sections. The first section deals with bistable multi vibrators in fixed bias and self biased configurations. The board is designed to study and experimentally check the most popular configurations of Transistor Multi vibrators circuits. The second section affords the study of monostable multi vibrators with fixed bias and self bias and the relevant controls. The third section is dedicated to the astable multi vibrators. The board is supplied complete with a set of stackable, plug •in cables of suitable lengths and colours and with a training manual. The various set, reset and trigger controls are analyzed. Power supply:  $\pm 15$  Vdc, 100 mA

#### **Technical Specification :**

##### **Experiments:**

Analysis of the monostable multi vibrator with fixed and self bias  
Analysis of the astable multi vibrator  
Check of the trigger control  
Analysis of the bistable multivibrator with fixed and self bias  
Check of set, reset and trigger controls operation.



**Laboratory Instrument India**

