



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Plunge And Fixed Base Variable Speed Router Kit

**Product Code :**  
LBNY-0004-770007



**Description :**

Plunge And Fixed Base Variable Speed Router Kit

**Technical Specification :**

Plunge And Fixed Base Variable Speed Router Kit

Electronic variable speed plunge and fixed base router kit; soft start, micro-fine depth adjustments  
2-1/4-horsepower router combination kit; variables speeds of 8,000 to 25,000 rpm Aluminum construction;  
wooden handles on fixed base; soft-grip handle on plunge base.

Includes plunge and fixed bases, collet chucks, wrenches, chip shields, carrying case 6.2 by 19-1/2 by 14-1/2 inches; 19-1/2 pounds

**Benchtop Router Table** Benchtop router table is compatible with most routers Includes 2 adjustable featherboards, dust collection port, dual outlet switch and 6-foot power cord Aluminum fence with MDF face plates and aluminum router mounting plate Includes mounting hardware for a variety of routers Work surface measures 27-inches by 18-inches

**Deluxe Router Edge Guide** Guides cut along the workpiece or up to 8" in from edge. Easily converts into circle guide for making arcs and circles up to 32" in diameter. Includes router guide, pivot plate, dust extraction hood, vacuum hose adaptor, attachment hardware.

---

# Laboratory Instrument India

