



Email : sales@laboratoryinstrumentindia.com

Product Name :
Chisel Mortiser Wood working Machine

Product Code :
LBNY-0004-770008



Description :

Chisel Mortiser Wood working Machine

Technical Specification :

Chisel Mortiser Wood working Machine
Workpiece width : 250mm
Mortising length : 200mm
Mortising depth : 110mm
Mortising width : 22mm
Spindle speed: 2800r/min
Motor power : 1.1kw
Overall dimension : 740x930x1350mm
Weight : 240kg



Laboratory Instrument India