



Email : sales@laboratoryinstrumentindia.com

Product Name :
Double Dowelling Machine

Product Code :
LBNY-0004-770009



Description :

Double Dowelling Machine

Technical Specification :

Double Dowelling Machine

Automatic depth control system allows for easy set-up for different material thicknesses.

Variable-spaced bushings allow a variety of positions for the pocket holes.

Self-adjusting clamp adjusts automatically to material thickness.

Repeatable clamp makes multiple pieces quick and easy .

Screw-length gauge recommends the appropriate screw length for material thickness.

Repeatable clamp makes multiple pieces quick and easy.

Built-In Drill Bit collar set up does not need to be changed for different material thickness.

Removable Dust Port can be hooked up to a Porter Cable Dust Extractor or standard vacuum.



Laboratory Instrument India